

EEGer⁴™

Neurofeedback Software

Technical Manual Version 4.3.0

This manual contains information intended for licensees of EEGer software.

All information is Copyright © 2013 EEG Software and all rights are reserved by EEG Software.

Initial release 31 July 2012

TM43001 16 August 2012 for standardized revision control

TM43002 09 December 2012 to add data acquisition methodology

TM43003 30 January 2013 for added explanation of amplifier test/data

TM43004 3 December 2013 for 430m release

Table of Contents

Computer Requirements.....	4
Timing.....	5
Filters.....	6
Filters Provided With EEGer4.....	6
Acquisition Components.....	10
Layout and Feedback Modes.....	13
Examples of modes.....	17
Data Storage Format.....	18
Appendix A: Filter Bandpass Characteristics.....	19
Appendix B: Device Bandpass Characteristics.....	31
Appendix C.....	39
Mapping of trace (numbers) to game strands.....	39
Table of modes by layout.....	46
Dataflow Diagrams.....	55
Appendix D: Data formats.....	215
Common Definitions:.....	215
Raw Data File Format.....	221
Summary File Format.....	227
Appendix E: Data Acquisition Methodology.....	234
Appendix F: Signal Generator.....	235
Technical Description.....	235
Description.....	236
Digital to Analog Converter (DAC).....	236
Attenuator.....	237
SigDriver.....	240
Frequency File Generator (FreqFileGen).....	243
Calibration process.....	245

Computer Requirements

EEGer software requires one or two computers to operate (depending upon user configuration selections). EEGer executes on the following operating system configurations:

- Windows XP 32-bit
- Windows Vista 32-bit
- Windows Vista 64-bit
- Windows 7 32-bit
- Windows 7 64-bit

Testing with pre-release Windows 8 showed EEGer performed correctly but final testing cannot be performed until actual release of Windows 8.

The most sensitive element in a computer system (for EEGer) is the graphics interface. Some graphics chipsets/drivers exhibit poor performance, causing apparent display lagging although acquisition and processing continue normally.

Recommended minimum computer requirements:

	Single computer system	Therapist computer	Client/Game computer
Processor speed	2 GHz	1.8 GHz	1.8 GHz
Memory	XP – 1 GB Vista/7 – 3 GB	XP – 512 KB Vista/7 – 2 GB	XP – 512 KB Vista/7 – 2 GB
Storage	250 GB	80 GB	80 GB
Video card/chipset At least DirectX 9.0c supported. Minimum resolution 1024x768.	Extended desktop support for an external monitor (and external monitor connector). High-level gaming performance. Note: ATI/AMD or nVidia recommended since not all Intel graphics have required performance.	512 MB memory with mid-level gaming performance	512 MB memory with mid-level gaming performance
Communication ports	USB for EEGer dongle+ USB/serial for acquisition device	USB for EEGer dongle+ USB/serial for acquisition device+ ethernet/serial for game connection link	ethernet/serial for therapist connection link

Timing

EEGer processes EEG samples at 256 Hz. Each “frame” (1/256 of a second) data is stored, filtered, and decision-tested. There are some inherent delays in the filtering process since multiple samples are needed to provide a filter “output”. The default timing/delays used in EEGer are as follows:

Sample acquisition timing

Although the nominal sample time is 4 milliseconds per sample, USB interfaces transmit data in blocks so there is a variable time based on block size. Range is 0 to about 32 milliseconds ‘lateness’ in each sample. This time disregards any acquisition component internal delays (settling times).

Filtering timing

This timing depends on the number of filter stages. EEGer default is 2 stages so the delay is 8 milliseconds for a signal to “exit” a filter. Actual computation time is less than a microsecond.

Decision logic

This is the time it takes for the software to smooth the raw cyclic data. It depends on user-selection of a smoothing value which ranges from 0.1 to 0.9 seconds. EEGer default is 0.5 seconds but actual delay depends on significance of new sample (larger signals have more impact on the smoothed value).

Clinician display

The clinician display process runs at a rate between 25 and 40 Hz so the data is displayed within 25 to 40 milliseconds of computation.

Transmission to client

For a single-computer system, data is transmitted using an internal TCP/IP transmission with a delay time of about 100 microseconds (until ready for receipt).

For a two-computer system with serial connections, data is sent at a 115,200 baud rate to the client computer. Each message ranges from 8 bytes to about 80 bytes so the maximum delay is about 10 milliseconds.

For a two-computer system with ethernet connections, the maximum delay is about a millisecond.

Client data processing

The client feedback/game software runs at a 40 millisecond interval. Aural cues are processed immediately. Visual data is smoothed at rates dependent on the display so that “strobe” effects are not generated. The maximum processing delay is thus 25 milliseconds.

Overall latency is thus the sum of delays ranging from each stage:

Best case: $0+8+1+25+\text{smoothing} = 34$ milliseconds plus any smoothing delays

Worst case: $32+8+10+25+\text{smoothing} = 75$ milliseconds plus any smoothing delays

Filters

Filters are characterized by many values. These typically include:

- a) Filter type (moving average, FIR, IIR, JTFA, wavelet, etc.)
- b) Number of stages (filter order)
- c) Rolloff (frequency) characteristics
- d) Measurement points (edge, corner, 50%, etc.)
- e) Ripple
- f) Impulse/step (transient) response
- g) Phase accuracy
- h) Delay

EEGer4 uses IIR digital filtering to separate out frequency bands of interest.

IIR (Infinite Impulse Response) filters are characterized by good amplitude fidelity but poor phase fidelity. FIR (Finite Impulse Response) filters are characterized by relatively poor amplitude fidelity but good phase fidelity. FIR filters typically require many more computations than IIR filters and consequently have a longer filter delay.

Also, there are many kinds of digital filter computations (Butterworth,), each with their own characteristics. The actual filter computation logic in EEGer4 is performed using biquad computations to model the required polynomials. Each “quad” contains the coefficients needed. Each quad corresponds to one filter order. This computation looks like this:

$$y_n = a0 * x_n + a1 * x_{n-1} + a2 * x_{n-2} - b1 * y_{n-1} - b2 * y_{n-2}$$

Equation 1

where x_n is the input sample and y_n is the output sample.

Filters Provided With EEGer4

There are three sets of filters provided with EEGer4, all sampled and computed at 256 Hz.

The first set consists of precomputed IIR elliptical filter coefficient sets in the range 0 to 50 Hz in 1/8 Hz steps. This set is the same set used in earlier versions of EEGer. The set was generated with specifications of:

ripple= 0.5db
rolloff= -60db (lowpass) or -30db (bandpass)
order= 1 (lowpass) or 2 (bandpass)

The second set of filters is similar (IIR elliptical filter coefficients) but with dynamically-computed values and a choice of order ranging from 2 to 5. Also, this filter set has a step size of 0.1 Hz if the high frequency is more than 1 Hz and 0.001 Hz if high frequency is less than 1 Hz.

ripple= 0.5db
rolloff= -30db
order= 2 to 8

The third set of filters is mechanized using an open-source filter module (FIDlib) which is also used by some other neurofeedback manufacturers. The specific filter types EEGer4 uses with FIDlib are BsBu (Bandstop Butterworth), BpBu (Bandpass Butterworth), and LpBu (Lowpass Butterworth). These filters all use the following specifications:

order= 2 to 8

width specified at -3db points (.707 of peak)

A comparison of the three methods is shown in Figures 1,2,3,4.




Figure 1




Figure 2




Figure 3

Here is a comparison of the filter sets for various orders:




Figure 4

Note that the EEGer4 frequency specifications are to the edges of the flat response while the FIDlib specifications are to the -3db points. Full bandpass characteristics are shown for each filter method are shown in Appendix A. Appendix E describes the methodology used to acquire this data.

Acquisition Components

EEGer4 supports amplifier/encoder components from many manufacturers. Each component has (or may have) different frequency response characteristics. The following encoder/amplifiers are supported by EEGer4 (Pass/Fail/Untested noted in the last column):

Device Name	Manufacturer	Hardware Interface	Interface method	Support status	EEG Channels Supported	P/F
BrainLynx	J&J Engineering	USB	Mfr DLL	Currently supported	2	P
C2mini	J&J Engineering	USB	Mfr DLL	Obsolete/ supported	2	U
C2	J&J Engineering	USB	Mfr DLL	Obsolete/ supported	2	U
C2+mini	J&J Engineering	USB	Mfr DLL	Currently supported	2	P
C2+	J&J Engineering	USB	Mfr DLL	Currently supported	2	P
Spectrum 2	J&J Engineering	USB	Mfr DLL	Currently supported	2	P
Spectrum 4	J&J Engineering	USB	Mfr DLL	Currently supported	2/4	P
esiPro 2.2	TeleDiagnostic Systems	USB	Mfr DLL	Currently supported	2	P
esiPro 4.3	TeleDiagnostic Systems	USB	Mfr DLL	Currently supported	2/4	P
A200	Phoenix Neuro Systems	USB	Mfr DLL	Currently supported	2	P
A400	Phoenix Neuro Systems	USB	Mfr DLL	Currently supported	2/4	P
ProComp2	Thought Technology	Serial port	Internal code (serial interface)	Currently supported	2	P
ProComp+	Thought Technology	Serial port	Internal code (serial interface)	Currently supported	2	P
Infiniti	Thought Technology	Serial port	Internal code (serial interface)	Currently supported	2	P
ProComp5	Thought Technology	Serial port	Internal code (serial interface)	Currently supported	2	U

Device Name	Manufacturer	Hardware Interface	Interface method	Support status	EEG Channels Supported	P/F
ProComp2	Thought Technology	USB	Mfr DLL	Currently supported	2	P
ProComp+	Thought Technology	USB	Mfr DLL	Currently supported	2/4	P
Infiniti	Thought Technology	USB	Mfr DLL	Currently supported	2/4	P
ProComp5	Thought Technology	USB	Mfr DLL	Currently supported	2/4	U
Brainmaster 2E	Brainmaster	Serial port	Mfr DLL	Currently supported	2	P
Atlantis 2	Brainmaster	USB with serial port	Mfr DLL	Currently supported	2	P
Atlantis 4	Brainmaster	USB with serial port	Mfr DLL	Untested/no device but same DLL as all Brainmaster devices	2/4	U
Discovery	Brainmaster	USB with serial port	Mfr DLL	Under test/no device	2/4 (24)	U
Pet2.0	Brainquiry	Bluetooth serial	Internal code (serial interface)	Obsolete/no device but previously supported	2	U
QPET	Brainquiry	Bluetooth serial	Mfr DLL	Obsolete/no device but previously supported	2	U
Pendant-EEG	Pocket-Neurobics	Bluetooth serial	Internal code (serial interface)	Under test/no device but previously supported	2	U
*-Wiz		USB/Bluetooth	Internal code	Supported	2/4	P

Device Name	Manufacturer	Hardware Interface	Interface method	Support status	EEG Channels Supported	P/F
(external amplifier)	-----	USB	A/D device supported by Measurement Computing Corporation InstaCal software	Currently supported (used for system testing)	2/4	P

The acquisition devices use a variety of interface methods, usually serial, Bluetooth, or USB connections. Each device manufacturer has a custom interface method. EEGer has interface modules written to support all the above devices using either direct programming or the device manufacturer's provided interface methodology (i.e. DLL or other interface program). All the devices from J&J Engineering use a common interface DLL, differing only by a command code, and a common EEGer interface module.. All the C2+ based devices (including BrainLynx and Spectrum 2/4) use a single command code (and a common printed circuit board). All the devices from Brainmaster use a common DLL and a single EEGer interface module. All the Thought Technology devices using serial interfaces use generic Windows serial interfaces and a single EEGer interface module. All the Thought Technology devices using the USB (TTUSB) interface use a common DLL and a single EEGer interface module. Since the interfaces are common between device models, extensive testing was only necessary for devices using each interface. A simpler test was used to confirm that the other models supported could also communicate across the interface.

Appendix B contains the bandpass characteristics of these components. These bandpass characteristics are in addition to the filter characteristics of whatever filter method and frequency band selected.

Appendix E contains the methodology used to acquire the test data for Appendix A and B.

Layout and Feedback Modes

EEGer4 provides a number of screen configurations (called “layouts”) of EEG data. Each layout has a number of “feedback modes” with predefined usage of each element of the screen. All layouts have two or four lowpass EEG traces at the top of the screen and some number of additonally filtered traces below.

The current list of layouts (and the short titles) is:

5= 5-trace
6i= 6-trace,inhibit
6r= 6-trace,reward
8= 8-trace
14= 14-trace
14m= 14-trace,monitor
7= 7-trace,2monitor
14a= ChanA,screen
14ab= 2-chan,screen
14b= ChanB,screen
5r= Reward-only
6m= 6-trace,monitor
7m= 7-trace,monitor
8z= 8-trace-4zcomp
8p= 8-trace-2-rewards
10c=10-trace,coherences
10r=10-trace,ratios

Each feedback mode (internally) consists of a sequence of filter operations and decisions tied to particular elements (traces both visible and invisible) of the layout. Appendix C has the block diagrams of each of these feedback modes. The operation blocks on those diagrams use the following operations/decisions (names shown in **bold**). Certain generic operations are used. Most importantly, smoothing is performed using an Exponentially-Weighted-Moving-Average (EWMA) filter characterized by

$$G_n = G_{n-1} * \frac{k_t - 1}{k_t} + g_n * \frac{1}{k_t}$$

Equation 2

Also, a threshold decision is made using a comparison between a short-term (EWMA) average and a user-specified threshold value.

lowpass

Lowpass filter and DC correct a raw input from an acquisition component.

The input value is compared to the current “artifact threshold” to determine if the data is a spike (short

term amplitude value) or artifact (longer term amplitude value). If the signal is determined to be an "artifact" (probably due to eyeblinks, muscle motion, etc.), a zero value is the output of the operation. Once in an artifact condition, that condition will be held for a short time after the condition disappears (allowing the data to stabilize). If the input is NOT an artifact, the input value is lowpass filtered. The output of the filter is DC-corrected (to center the data display) and is the output of this operation. It is also integrated using an EWMA filter with a short-term, user-specified time constant (typically 0.5 seconds, with a range of 0.1 to 0.9 seconds).

bandpass

Bandpass filter a lowpass data sample and check for target threshold.

The input value is bandpass filtered. The signal is integrated using an EWMA for comparison with the target threshold. This block has several outputs: the signal, the 'average' value (actually the effective Peak-Peak voltage), and an integral-value-over-threshold flag. There are actually four submodes of this (fundamental) operation depending on the elements/traces configured:

OneChannelReward – if element is a reward element and only one input was specified
 SumReward --if element is a reward element and two inputs were specified
 (the two input values are arithmetically added together before filtering, etc.)
 OneChannelInhibit -- if element is an inhibit element and only one input was specified
 TwoChannellInhibit ----if element is an inhibit element and two inputs were specified
 (the two input values are arithmetically added together before filtering, etc.)

differ

Subtract the second specified element from the first specified element, then proceed as in OneChannelReward bandpass.

psync

Coherence measure between peak values.

The two (specified) input streams are each narrow-band filtered (using the reward frequency band limits) yielding two values (x and y). A cross correlation is then performed on a window (w) of the value histories:

$$\frac{\int X(t) Y(t)}{\sqrt{\int X(t) X(t) \int Y(t) Y(t)}}$$

Equation 3

This reduces in practice to

$$\frac{\sum_1^W X(t) Y(t)}{\sqrt{\sum_1^W X(t) X(t) * \sum_1^W Y(t) Y(t)}}$$

Equation 4

The window width is user-specifiable but defaults to 0.5 seconds. The resulting correlation value is smoothed using the standard EWMA filter. This value ranges from 0 to 1. It is multiplied by a constant (user specified but default to 10.0) to place it in an appropriate display range (0-10) for typical EEG displays. The scaled value is the principal output of this operation.

async

Coherence measure between wave slope angles.

The two (specified) input streams are each narrow-band filtered (using the reward frequency band limits) yielding two values (x and y). Each stream history is examined (looking 'backward' for one cycle) to find the min/max of the wave. The difference between the (smoothed) average amplitude of the stream and the current value is computed. A cross correlation is performed of the difference values using Equation 4. The output of this is the async value which ranges from 0 to 2. It is multiplied by a constant (user specified but default to 10.0) to place it in an appropriate display range (0-10) for typical EEG displays. The scaled value is the principal output of this operation.

gasync

Coherence measure exactly like async except using the lowpass (wideband) data instead of the narrow band data.

aminusb

Compute difference between short-term average of two channels.

This operation filters each channel of data separately and computes a short-term moving average value. The two short-term values are subtracted and the result added to the threshold value to give an output (display) value.

pdelta

Compute variation in timing between peaks on two channels.

This routine determines the 'time' of the peak value of the last cycle's samples in each channels. The smoothed time is compared to the smoothed time of the other channels and a normalized delta time value is determined. The delta time is smoothed and compared to the instantaneous delta time and used as the output value. The value is subtracted from 1 so that zero variation in peak times results in a perfect correlation factor. It is multiplied by a constant (user specified but default to 10.0) to place it in an appropriate display range (0-10) for typical EEG displays. The scaled value is the principal output of this operation.

ratio

Compute the ratio of two channels or two streams of data.

Depending in the input specifications, this operation determine the ratio of the inputs. It is a user configuration selection whether the ratio is squared (for power) or not (for voltage).

diffsum

This operation computes the difference of the filtered samples divided by the sum.

$$\frac{A_n - B_n}{A_n + B_n}$$

Equation 5

unity

This operation is very similar to the diffsum operation except the value is subtracted from 1.

$$1 - \frac{A_n - B_n}{A_n + B_n}$$

Equation 6

zcomp

This operation is not a filtering operation but examines various output combinations of data computed by the Applied NeuroScience Inc. (ANI) zscore module. The zcomp logic examines each user-specified element of the ANI data outputs and tracks the percentage of rewardable states. The composite percentage is the output of this operation block. Further details of the comparison logic are described in the Operator Manual.

qavg

This operation computes the average Psync values between all the pairs of channels with respect to one channel.

qpsfun

This operation has several sub functions. Basically, for each of the currently-matching filter streams (each with matching frequency limits and operating in the A_Psync feedback mode), the (A_Psync) values are saved. The resulting operation depends on the submode.

QPSavg:

The values are summed and divided by the number of streams.

QPSlag:

The values are integrated using an EWMA filter.

QPSdev:

The deviation of the values is computed.

The following operations are not filtering operations but reward decision operations that receive inhibit and reward inputs and make various decisions about scoring, sounds, etc.

bsmr

This operation block accepts the inhibit and reward threshold comparisons and determines rewardable state for SMR and EXP protocol classes. It also determines if a reward event is to occur based on time-on-task, event rate, time between rewards, and other controlling criteria.

at

This operation block accepts the inhibit and reward threshold comparisons and determines rewardable state for the AT protocol class. Reward sounds are determined in the feedback games/displays. The feedback consists of two environmental tones signaling the general *range* of the client's state-of-relaxation, and tones that signal momentary bursts of EEG in one or the other reward band. The alpha environment is a running stream, and its tone is a high-pitched bell. The theta environment is ocean waves, and the theta tone is a low-pitched bell. As relaxation deepens, the ongoing environmental sound will cross-fade from stream to ocean and low-pitched bell tones will replace the high-pitched tones. In cases where both frequency ranges are over threshold at the same time, the system is biased toward signaling the alpha state. The user has control over the "fade" rate when transitioning between states. If neither signal is over threshold, there will be no change in relative volumes (fading).

The following operations are not strictly operations performed by EEGer4. They consist of various selection of outputs from the Applied NeuroScience Inc. (ANI) zscore module for display/reward controls.

zasymm - asymmetry measure

zcohere - coherence measure

zphase - absolute phase measure

zabspa - absolute power

zabspb - absolute power

zrelpa - relative power

zrelpb - relative power

zpratioa – power ratios

zpratiob – power ratios

For an explanation of these computations, please see the applicable ANI documentation.

All the feedback modes are diagrammed in Appendix C.

Examples of modes

Some explanation of (internal) entries in the samples below:

display means display the output

game specifies the game strand

extra enables multiple reward modes

proc gives text displayed to therapist

kind is the layout code where C means channel, I means Inhibit, R means reward.

Each stream has its own frequency bandwidth for bandpass filtering.

A simple example of the mode logic (all diagrammed in Appendix C) is Differ.

kind=CCIRI

10=lowpass,in=0,out=0,display

>Lowpass channel A, display
on trace 0

20=lowpass,in=1,out=1,display

>Lowpass channel B, display
on trace 1

30=bandpass,in=0-1,out=2,game=0,display

>Bandpass trace 0 and 1
added, display on trace 2

40=bandpass,in=0-1,out=4,game=3,display

>Bandpass trace 0 and 1

50=differ,in=0-1,out=3,game=1,display,extra,proc=Diff	added, display on trace 4 >Subtract trace 1 from trace 0, bandpass, display on trace 3
60=bsmr,in=3-3,inhibit=2-4	>Reward if trace 3 rewardable and no inhibit on traces 2 and 4

Another example is the 8-trace Dual mode:



kind=CCIRIIIRI	
10=lowpass,in=0,out=0,display	
20=lowpass,in=1,out=1,display	
30=bandpass,in=0,out=2,game=0,display	
40=bandpass,in=0,out=4,game=3,display	
50=bandpass,in=0,out=3,game=1,display,extra,proc=Ampl	> note inputs come from a single channel
51=bandpass,in=1,out=5,game=4,display	
52=bandpass,in=1,out=7,game=7,display	
53=bandpass,in=1,out=6,game=6,display,extra,proc=Ampl	> note inputs come from a single channel
60=bsmr,in=3-6,inhibit=2-4-5-7	> both rewards and all 4 inhibits participate


Data Storage Format



The internal data structures for raw and summary data are described by the header files listed in Appendix D. Further information and guidance can be obtained from EEG Software on request.



Appendix A: Filter Bandpass Characteristics



All these data results are based on a 40 microvolt peak-peak input voltage.



































All lowpass filtering (0 to 30,40,50 Hz) use the same filters as Dynamic 2 settings.






















Appendix B: Device Bandpass Characteristics





Bandpass characteristics plots based on 40 microvolt peak-peak inputs:




















Other graphs and plots available from EEG Software as requested.

Device	Serial number	Notes
Spectrum 2	12399	
Spectrum 4	-----	Prototype without serial number
Brainmaster 2E	4290	
Atlantis II	40403	
esiPro 2.2	2C099132	
esiPro 4.3	B117014	
ProComp2	BC1063	
ProComp+	AD3526	Flexpro=G3936
Infiniti		
Pendant 9v		
Pet2.0		
QPET		
BrainLynx		

Appendix C

Mapping of trace (numbers) to game strands

These tables show the connection between the trace data (left column) and the game strand that the data is sent to.

This data is primarily of use for the sound specialization and for developers of new games.

(5) 5-trace 2 chan,1 reward,2 inhibit

2	0
4	3
3	1

(6i) 6-trace,inhibit 2 chan,1 reward,3 inhibit

2	0
4	3
3	1
5	7

(6r) 6-trace,reward 2 chan,2 reward,2 inhibit

2	0
5	3
3	1
4	5

(8) 8-trace 2 chan,2 reward,4 inhibit

2	0
4	3
3	1
5	4
7	7
6	5

(10c) 10-trace,coherences 4 chan,2 reward,4 inhibit

4	0
6	3
5	1
7	4
9	7
8	5

(10r) 10-trace,ratios 4 chan,2 reward,4 inhibit

4	0
6	3
5	1
7	4
9	7
8	5

(12) 12-trace 4 chan,4 monitor,1 reward,3 inhibit

8	0
9	4
11	3
10	1

(14a) ChanA,screen 14 traces,no rewards

(14ab) 2-chan,screen 14 traces,no rewards

(14b) ChanB,screen 14 traces,no rewards

(5r) Reward-only 2 chan,1 reward,2 monitor

4	1
---	---

(8p) 8-trace 4 chan,2 reward,2 inhibit

4	0
5	1
6	5
7	3

(5) 5-trace 2 chan,2 reward,1 inhibit AlphaTheta

2	0
4	3
3	1

(6i) 6-trace,inhibit 2 chan,1 reward,3 inhibit AlphaTheta

2	0
4	3
3	1
5	7

(8) 8-trace 2 chan,2 reward,4 inhibit AlphaTheta

2	0
4	3
3	1
5	7

(8q) 8-trace 4 chan,1 reward,3 inhibit

4	0
5	4
7	3
6	1

EEGer4 Technical Manual

(8z) 8-trace-4zcomp 4 chan,4 reward

4	0
5	4
6	1
7	3

(10c) 10-trace,coherences 4 chan,2 reward,4 inhibit

4	0
6	3
5	1
7	4
9	7
8	5

(10r) 10-trace,ratios 4 chan,2 reward,4 inhibit

4	0
6	3
5	1
7	4
9	7
8	5

(12) 12-trace 4 chan,4 monitor,1 reward,3 inhibit

8	0
9	4
11	3
10	1

Alphabetical Index

ABDIFFABDIFF (DIFFAB + DIFFAB)	79
ABSUMABDIFF (SUMAB + DIFFAB)	80
ABSUMABSUM (SUMAB + SUMAB)	78
AminusB (A channel relationship to B channel)	62, 88, 144
Async (comodulation measure between channel A and B)	60, 86, 142
BminusA (B channel relationship to A channel)	63, 89, 145
D/S-Ratio (Ratio of A-B to A+B)	95
D/S-Ratio (Ratio of AB diff/AB sum)	69
DAsync (4 channel AsyncAB+AsyncCD)	181
DAsyncBPAB (4 channel AsyncAB+BP C)	179
DDiffABDiffCD (4 channel DiffAB DiffCD)	173
DDiffABSumCD (4 channel DiffAB SumCD)	172
DDiffSum (4 channel Diff-SumAB+Diff-SumCD)	199
DDiffSumBPAB (4 channel Diff-SumAB+BP C)	197
DDiffSumBPCD (4 channel Diff-SumCD+BP A)	198
DGAsync (4 channel GAsyncAB+GAsyncCD)	184
DGAsyncBPAB (4 channel GAsyncAB+BP C)	182
DGAsyncBPCD (4 channel GAsyncCD+BP A)	183
Differ (channel A minus channel B)	58, 84, 101, 128, 132, 136, 140
DMinus (4 channel MinusAB+MinusCD)	190
DMinusBPAB (4 channel MinusAB+BP C)	188
DMinusBPCD (4 channel MinusCD+BP A)	189
DPdelta (4 channel PsyncAB+PsyncCD)	187
DPdeltaBPAB (4 channel PsyncAB+BP C)	185
DPdeltaBPCD (4 channel PsyncCD+BP A)	186
DPsync (4 channel PsyncAB+PsyncCD)	178
DPsyncABCD (4 channel two psync)	167
DPsyncAvg2 (4 channel 2 averaged psyncs)	169
DPsyncAvgBP (4 channel averaged psync + A)	168

EEGer4 Technical Manual

DPsyncBPAB (4 channel PsyncAB+BP C)	176
DPsyncBPCD (4 channel PsyncCD+BP A)	177
DRatio (4 channel RatioAB+RatioCD)	196
DRatioBPAB (4 channel RatioAB+BP C)	194
DRatioBPCD (4 channel RatioCD+BP A)	195
DSingleA (4 channel single A)	163
DSingleAB (4 channel single A B)	174
DSingleB (4 channel single B)	164
DSingleC (4 channel single C)	165
DSingleCD (4 channel single C D)	175
DSingleD (4 channel single D)	166
DSsyncBPCD (4 channel AsyncCD+BP A)	180
DSumABDiffCD (4 channel SumAB DiffCD)	171
DSumABSumCD (4 channel SumAB SumCD)	170
Dual (two channels of data input)	72
DUnity (4 channel UnityAB+UnityCD)	202
DUnityBPAB (4 channel UnityAB+BP C)	200
DUnityBPCD (4 channel UnityCD+BP A)	201
DZcomp (4 channel ZcompAB+ZcompCD)	193
DZcomp (4 channel zcomposite)	160
DZcompBPAB (4 channel ZcompAB+BP C)	191
DZcompBPAB (4 channel zcompositeAB + single A)	161
DZcompBPCD (4 channel ZcompCD+BP A)	192
DZcompBPCD (4 channel zcompositeCD + single A)	162
GAsync (global comodulation measure between channel A and B)	61, 87, 143
MonitorA (1-ch + 12 monitors)	103
MonitorAB (2-ch + 12 monitors)	104
MonitorB (1-ch + 12 monitors)	105
PDelta (B channel relationship to A channel)	90
PDelta (Peak Time Coherence)	64, 146
Psync (synchrony measure between channel A and B)	59, 75, 85, 141
PsyncABAB (comodulation measure AB twice)	102

EEGer4 Technical Manual

QABCD (4 channel amplitude)	158
QPSAvg (4 channel Average Psyncs)	159
QPSAvg (4 channel PSync Averaged Reward + 3 Inhibits)	208
QPSAvgA (4 channel PSync Averaged + 3 Inhibits)	203
QPSdev (4 channel PSync Deviation Reward + 3 Inhibits)	210
QPSlag (4 channel PSync Smoothed Reward + 3 Inhibits)	209
QQSingleB (4 channel Monitor + Reward+Inhibits)	211p.
QQSingleC (4 channel Monitor + Reward+Inhibits)	213
QQSingleD (4 channel Monitor + Reward+Inhibits)	214
QSingleA (4 channel BP-A)	204
QSingleB (4 channel BP-B)	205
QSingleC (4 channel BP-C)	206
QSingleD (4 channel BP-D)	207
QZcomp (4 channel zcomposite)	157
RatioA (Ratio of 2 streams from 1 channel)	91
RatioA (Ratio of rwd/inhib)	65
RatioAB (Ratio of A to B in reward band)	67, 93
RatioB (Ratio of 2 streams from 1 channel)	92
RatioB (Ratio of rwd/inhib)	66
RatioBA (Ratio of B to A in reward band)	68, 94
SingleA (one channel of data input)	55, 81, 98, 125, 129, 133, 137
SingleA (SingleA)	73
SingleB (one channel of data input)	56, 82, 99, 126, 130, 134, 138
SingleB (SingleB)	74
Sum (sum of two channels of data input)	57, 83, 100, 127, 131, 135, 139
Unity (1-Ratio of A-B to A+B)	96
Unity (1-Ratio of AB diff/AB sum)	70
ZAbsPwrA (Zscore Abs Amp A)	110, 119, 151
ZAbsPwrB (Zscore Abs Amp B)	111, 120, 152
ZAsymm (Zscore amplitude asymmetry)	107, 116, 148
ZCoher (Zscore coherence)	108, 117, 149
ZCompAB (Zcomposite + Zcomposite)	77

EEGer4 Technical Manual

ZCompAB (ZComposite)	71, 97, 106, 147
ZCompBPAB (Zcomposite + bandpass)	76
ZPhase (Zscore phase)	109, 118, 150
ZPRatioA (Zscore Power ratio A)	114, 123, 155
ZPRatioB (Zscore Power ratio B)	115, 124, 156
ZRelPwrA (Zscore Rel Pwr A)	112, 121, 153
ZRelPwrB (Zscore Rel Pwr B)	113, 122, 154

Table of modes by layout

SMR/EXP 5-trace

100 SingleA (one channel of data input)	CCIRI (5)	55
101 SingleB (one channel of data input)	CCIRI (5)	56
110 Sum (sum of two channels of data input)	CCIRI (5)	57
111 Differ (channel A minus channel B)	CCIRI (5)	58
120 Psync (synchrony measure between channel A and B)	CCIRI (5)	59
130 Async (comodulation measure between channel A and B)	CCIRI (5)	60
140 GAsync (global comodulation measure between channel A and B)	CCIRI (5)	61
150 AminusB (A channel relationship to B channel)	CCIRI (5)	62
151 BminusA (B channel relationship to A channel)	CCIRI (5)	63
152 PDelta (Peak Time Coherence)	CCIRI (5)	64
153 RatioA (Ratio of rwd/inhib)	CCIRI (5)	65
154 RatioB (Ratio of rwd/inhib)	CCIRI (5)	66
155 RatioAB (Ratio of A to B in reward band)	CCIRI (5)	67
156 RatioBA (Ratio of B to A in reward band)	CCIRI (5)	68
158 D/S-Ratio (Ratio of AB diff/AB sum)	CCIRI (5)	69
159 Unity (1-Ratio of AB diff/AB sum)	CCIRI (5)	70
180 ZCompAB (ZComposite)	CCIRI (5)	71
700 ZAsymm (Zscore amplitude asymmetry)	CCIRI (5)	107
701 ZCohere (Zscore coherence)	CCIRI (5)	108

EEGer4 Technical Manual

702 ZPhase (Zscore phase)	CCIRI (5)	109
703 ZAbsPwrA (Zscore Abs Amp A)	CCIRI (5)	110
704 ZAbsPwrB (Zscore Abs Amp B)	CCIRI (5)	111
705 ZRelPwrA (Zscore Rel Pwr A)	CCIRI (5)	112
706 ZRelPwrB (Zscore Rel Pwr B)	CCIRI (5)	113
707 ZPRatioA (Zscore Power ratio A)	CCIRI (5)	114
708 ZPRatioB (Zscore Power ratio B)	CCIRI (5)	115

SMR/EXP 8-trace

200 Dual (two channels of data input)	CCIRIIRI (8)	72
210 SingleA (SingleA)	CCIRIIRI (8)	73
211 SingleB (SingleB)	CCIRIIRI (8)	74
221 Psync (synchrony measure between channel A and B)	CCIRIIRI (8)	75
260 ZCompBPAB (Zcomposite + bandpass)	CCIRIIRI (8)	76
261 ZCompAB (Zcomposite + Zcomposite)	CCIRIIRI (8)	77
270 ABSUMABSUM (SUMAB + SUMAB)	CCIRIIRI (8)	78
271 ABDIFFABDIFF (DIFFAB + DIFFAB)	CCIRIIRI (8)	79
272 ABSUMABDIFF (SUMAB + DIFFAB)	CCIRIIRI (8)	80

SMR/EXP 6,inhibit

400 SingleA (one channel of data input)	CCIIRI (6i)	81
401 SingleB (one channel of data input)	CCIIRI (6i)	82
410 Sum (sum of two channels of data input)	CCIIRI (6i)	83

EEGeri4 Technical Manual

411 Differ (channel A minus channel B)	CCIIIRI (6i)	84
420 Psync (synchrony measure between channel A and B)	CCIIIRI (6i)	85
430 Async (comodulation measure between channel A and B)	CCIIIRI (6i)	86
440 GAsync (global comodulation measure between channel A and B)	CCIIIRI (6i)	87
450 AminusB (A channel relationship to B channel)	CCIIIRI (6i)	88
451 BminusA (B channel relationship to A channel)	CCIIIRI (6i)	89
452 PDelta (B channel relationship to A channel)	CCIIIRI (6i)	90
453 RatioA (Ratio of 2 streams from 1 channel)	CCIIIRI (6i)	91
454 RatioB (Ratio of 2 streams from 1 channel)	CCIIIRI (6i)	92
455 RatioAB (Ratio of A to B in reward band)	CCIIIRI (6i)	93
456 RatioBA (Ratio of B to A in reward band)	CCIIIRI (6i)	94
458 D/S-Ratio (Ratio of A-B to A+B)	CCIIIRI (6i)	95
459 Unity (1-Ratio of A-B to A+B)	CCIIIRI (6i)	96
480 ZCompAB (ZComposite)	CCIIIRI (6i)	97
710 ZAsymm (Zscore amplitude asymmetry)	CCIIIRI (6i)	116
711 ZCohere (Zscore coherence)	CCIIIRI (6i)	117
712 ZPhase (Zscore phase)	CCIIIRI (6i)	118
713 ZAbsPwrA (Zscore Abs Amp A)	CCIIIRI (6i)	119
714 ZAbsPwrB (Zscore Abs Amp B)	CCIIIRI (6i)	120
715 ZRelPwrA (Zscore Rel Pwr A)	CCIIIRI (6i)	121
716 ZRelPwrB (Zscore Rel Pwr B)	CCIIIRI (6i)	122
717 ZPRatioA (Zscore Power ratio A)	CCIIIRI (6i)	123

EEGer4 Technical Manual

718 ZPRatioB (Zscore Power ratio B)	CCIIIRI (6i)	124
-------------------------------------	--------------	-----

SMR,EXP 6,reward

602 SingleA (one channel of data input)	CCIRRI (6r)	98
603 SingleB (one channel of data input)	CCIRRI (6r)	99
612 Sum (sum of two channels of data input)	CCIRRI (6r)	100
613 Differ (channel A minus channel B)	CCIRRI (6r)	101
634 PsyncABAB (comodulation measure AB twice)	CCIRRI (6r)	102
680 ZCompAB (ZComposite)	CCIRRI (6r)	106

EXP monitors

668 MonitorA (1-ch + 12 monitors)	CCMMMMMMMMMMMM (14a)	103
669 MonitorAB (2-ch + 12 monitors)	CCMMMMMMMMMMMM (14ab)	104
670 MonitorB (1-ch + 12 monitors)	CCMMMMMMMMMMMM (14b)	105

AT 5-trace

1100 SingleA (one channel of data input)	CCRRI (5)	125
1101 SingleB (one channel of data input)	CCRRI (5)	126
1110 Sum (sum of two channels of data input)	CCRRI (5)	127
1111 Differ (channel A minus channel B)	CCRRI (5)	128

AT 8-trace

1200 SingleA (one channel of data input)	CCRRIIRI (8)	129
1201 SingleB (one channel of data input)	CCRRIIRI (8)	130
1210 Sum (sum of two channels of data input)	CCRRIIRI (8)	131

EEGer4 Technical Manual

1211 Differ (channel A minus channel B)	CCRRIIRI (8)	132
---	--------------	-----

AT 6,inhibit

1300 SingleA (one channel of data input)	CCRRII (6i)	133
1301 SingleB (one channel of data input)	CCRRII (6i)	134
1310 Sum (sum of two channels of data input)	CCRRII (6i)	135
1311 Differ (channel A minus channel B)	CCRRII (6i)	136

SMR,EXP 5-trace reward-only

2000 SingleA (one channel of data input)	CCMMR (5r)	137
2001 SingleB (one channel of data input)	CCMMR (5r)	138
2010 Sum (sum of two channels of data input)	CCMMR (5r)	139
2011 Differ (channel A minus channel B)	CCMMR (5r)	140
2020 Psync (synchrony measure between channel A and B)	CCMMR (5r)	141
2030 Async (comodulation measure between channel A and B)	CCMMR (5r)	142
2040 GAsync (global comodulation measure between channel A and B)	CCMMR (5r)	143
2050 AminusB (A channel relationship to B channel)	CCMMR (5r)	144
2051 BminusA (B channel relationship to A channel)	CCMMR (5r)	145
2052 PDelta (Peak Time Coherence)	CCMMR (5r)	146
2060 ZCompAB (ZComposite)	CCMMR (5r)	147
2100 ZAsymm (Zscore amplitude asymmetry)	CCMMR (5r)	148
2101 ZCohere (Zscore coherence)	CCMMR (5r)	149
2102 ZPhase (Zscore phase)	CCMMR (5r)	150

EEGer4 Technical Manual

2103 ZAbsPwrA (Zscore Abs Amp A)	CCMMR (5r)	151
2104 ZAbsPwrB (Zscore Abs Amp B)	CCMMR (5r)	152
2105 ZRelPwrA (Zscore Rel Pwr A)	CCMMR (5r)	153
2106 ZRelPwrB (Zscore Rel Pwr B)	CCMMR (5r)	154
2107 ZPRatioA (Zscore Power ratio A)	CCMMR (5r)	155
2108 ZPRatioB (Zscore Power ratio B)	CCMMR (5r)	156

6000 QZcomp (4 channel zcomposite)	CCCCRRRR (8z)	157
6010 QABCD (4 channel amplitude)	CCCCRRRR (8z)	158
6020 QPSAvg (4 channel Average Psyncs)	CCCCRRRR (8z)	159

6100 DZcomp (4 channel zcomposite)	CCCCIRRI (8p)	160
6110 DZcompBPAB (4 channel zcompositeAB + single A)	CCCCIRRI (8p)	161
6114 DZcompBPCD (4 channel zcompositeCD + single A)	CCCCIRRI (8p)	162
6120 DSsingleA (4 channel single A)	CCCCIRRI (8p)	163
6130 DSsingleB (4 channel single B)	CCCCIRRI (8p)	164
6140 DSsingleC (4 channel single C)	CCCCIRRI (8p)	165
6150 DSsingleD (4 channel single D)	CCCCIRRI (8p)	166
6160 DPsyncABCD (4 channel two psync)	CCCCIRRI (8p)	167
6170 DPsyncAvgBP (4 channel averaged psync + A)	CCCCIRRI (8p)	168
6180 DPsyncAvg2 (4 channel 2 averaged psyncs)	CCCCIRRI (8p)	169

EEGer4 Technical Manual

6182 DSumABSumCD (4 channel SumAB SumCD)	CCCCIRRI (8p)	170
6184 DSumABDiffCD (4 channel SumAB DiffCD)	CCCCIRRI (8p)	171
6186 DDifftABSumCD (4 channel DiffAB SumCD)	CCCCIRRI (8p)	172
6188 DDifftABDiffCD (4 channel DiffAB DiffCD)	CCCCIRRI (8p)	173

10-trace,coherences

6200 DSsingleAB (4 channel single A B)	CCCCIRIIRI (10c,10r)	174
6210 DSsingleCD (4 channel single C D)	CCCCIRIIRI (10c,10r)	175
6370 DZcompBPAB (4 channel ZcompAB+BP C)	CCCCIRIIRI (10c,10r)	191
6380 DZcompBPCD (4 channel ZcompCD+BP A)	CCCCIRIIRI (10c,10r)	192
6390 DZcomp (4 channel ZcompAB+ZcompCD)	CCCCIRIIRI (10c,10r)	193
6220 DPsyncBPAB (4 channel PsyncAB+BP C)	CCCCIRIIRI (10c)	176
6230 DPsyncBPCD (4 channel PsyncCD+BP A)	CCCCIRIIRI (10c)	177
6240 DPsync (4 channel PsyncAB+PsyncCD)	CCCCIRIIRI (10c)	178
6250 DAsyncBPAB (4 channel AsyncAB+BP C)	CCCCIRIIRI (10c)	179
6260 DSsyncBPCD (4 channel AsyncCD+BP A)	CCCCIRIIRI (10c)	180
6270 DAsync (4 channel AsyncAB+AsyncCD)	CCCCIRIIRI (10c)	181
6280 DGAsyncBPAB (4 channel GAsyncAB+BP C)	CCCCIRIIRI (10c)	182
6290 DGAsyncBPCD (4 channel GAsyncCD+BP A)	CCCCIRIIRI (10c)	183
6300 DGAsync (4 channel GAsyncAB+GAsyncCD)	CCCCIRIIRI (10c)	184
6310 DPdeltaBPAB (4 channel PsyncAB+BP C)	CCCCIRIIRI (10c)	185
6320 DPdeltaBPCD (4 channel PsyncCD+BP A)	CCCCIRIIRI (10c)	186
6330 DPdelta (4 channel PsyncAB+PsyncCD)	CCCCIRIIRI (10c)	187

EEGer4 Technical Manual

6340 DMinusBPAB (4 channel MinusAB+BP C)	CCCCIRIIRI (10c)	188
6350 DMinusBPCD (4 channel MinusCD+BP A)	CCCCIRIIRI (10c)	189
6360 DMinus (4 channel MinusAB+MinusCD)	CCCCIRIIRI (10c)	190

10-trace, ratios


6200 DSsingleAB (4 channel single A B)	CCCCIRIIRI (10c,10r)	174
6210 DSsingleCD (4 channel single C D)	CCCCIRIIRI (10c,10r)	175
6370 DZcompBPAB (4 channel ZcompAB+BP C)	CCCCIRIIRI (10c,10r)	191
6380 DZcompBPCD (4 channel ZcompCD+BP A)	CCCCIRIIRI (10c,10r)	192
6390 DZcomp (4 channel ZcompAB+ZcompCD)	CCCCIRIIRI (10c,10r)	193
6420 DRatioBPAB (4 channel RatioAB+BP C)	CCCCIRIIRI (10r)	194
6430 DRatioBPCD (4 channel RatioCD+BP A)	CCCCIRIIRI (10r)	195
6440 DRatio (4 channel RatioAB+RatioCD)	CCCCIRIIRI (10r)	196
6450 DDifSumBPAB (4 channel Diff-SumAB+BP C)	CCCCIRIIRI (10r)	197
6460 DDifSumBPCD (4 channel Diff-SumCD+BP A)	CCCCIRIIRI (10r)	198
6470 DDifSum (4 channel Diff-SumAB+Diff-SumCD)	CCCCIRIIRI (10r)	199
6480 DUnityBPAB (4 channel UnityAB+BP C)	CCCCIRIIRI (10r)	200
6490 DUnityBPCD (4 channel UnityCD+BP A)	CCCCIRIIRI (10r)	201
6500 DUnity (4 channel UnityAB+UnityCD)	CCCCIRIIRI (10r)	202

6600 QPSAvgA (4 channel PSync Averaged + 3 Inhibits)	CCCCIRI (8q)	203
6610 QSingleA (4 channel BP-A)	CCCCIRI (8q)	204

EEGer4 Technical Manual


6620 QSingleB (4 channel BP-B)	CCCCIIRI (8q)	205
6630 QSingleC (4 channel BP-C)	CCCCIIRI (8q)	206
6640 QSingleD (4 channel BP-D)	CCCCIIRI (8q)	207

6800 QPSAvg (4 channel PSync Averaged Reward + 3 Inhibits)	CCCCMMMMIIRI (12)	208
6810 QPSlag (4 channel PSync Smoothed Reward + 3 Inhibits)	CCCCMMMMIIRI (12)	209
6820 QPSdev (4 channel PSync Deviation Reward + 3 Inhibits)	CCCCMMMMIIRI (12)	210
6840 QQSsingleB (4 channel Monitor + Reward+Inhibits)	CCCCMMMMIIRI (12)	211
6842 QQSsingleB (4 channel Monitor + Reward+Inhibits)	CCCCMMMMIIRI (12)	212
6844 QQSsingleC (4 channel Monitor + Reward+Inhibits)	CCCCMMMMIIRI (12)	213
6846 QQSsingleD (4 channel Monitor + Reward+Inhibits)	CCCCMMMMIIRI (12)	214

Dataflow Diagrams


100 SingleA (one channel of data input) CCIRI (5)

SingleA (one channel of data input)




101 SingleB (one channel of data input) CCIRI (5)

SingleB (one channel of data input)




110 Sum (sum of two channels of data input) CCIRI (5)

Sum (sum of two channels of data input)




111 Differ (channel A minus channel B) CCIRI (5)

Differ (channel A minus channel B)




120 Psync (synchrony measure between channel A and B) CCIRI (5)

Psync (synchrony measure between channel A and B)




130 Async (comodulation measure between channel A and B) CCIRI (5)

Async (comodulation measure between channel A and B)




140 GAsync (global comodulation measure between channel A and B) CCIRI (5)

GAsync (global comodulation measure between channel A and B)




150 AminusB (A channel relationship to B channel) CCIRI (5)

AminusB (A channel relationship to B channel)




151 BminusA (B channel relationship to A channel) CCIRI (5)


BminusA (B channel relationship to A channel)



PDelta (Peak Time Coherence)




RatioA (Ratio of rwd/inhib)




154 RatioB (Ratio of rwd/inhib) CCIRI (5)

RatioB (Ratio of rwd/inhib)




155 RatioAB (Ratio of A to B in reward band) CCIRI (5)

RatioAB (Ratio of A to B in reward band)




156 RatioBA (Ratio of B to A in reward band) CCIRI (5)

RatioBA (Ratio of B to A in reward band)




158 D/S-Ratio (Ratio of AB diff/AB sum) CCIRI (5)

D/S-Ratio (Ratio of AB diff/AB sum)



159 Unity (1-Ratio of AB diff/AB sum) CCI RI (5)


Unity (1-Ratio of AB diff/AB sum)



180 ZCompAB (ZComposite) CCIRI (5)


ZCompAB (ZComposite)

EEGer4 Technical Manual




Dual (two channels of data input)

EEGer4 Technical Manual





SingleA (SingleA)



SingleB (SingleB)


EEGer4 Technical Manual






ZCompBPAB (Zcomposite + bandpass)

EEGer4 Technical Manual




ZCompAB (Zcomposite + Zcomposite)

EEGer4 Technical Manual




EEGer4 Technical Manual



271 ABDIFFABDIFF (DIFFAB + DIFFAB) CCIRIIRI (8)


ABDIFFABDIFF (DIFFAB + DIFFAB)


EEGer4 Technical Manual




ABSUMABDIFF (SUMAB + DIFFAB)

EEGer4 Technical Manual









410 Sum (sum of two channels of data input) CCIIRI (6i)

Sum (sum of two channels of data input)




411 Differ (channel A minus channel B) CCIIRI (6i)

Differ (channel A minus channel B)




420 Psync (synchrony measure between channel A and B) CCIIRI (6i)

Psync (synchrony measure between channel A and B)




430 Async (comodulation measure between channel A and B) CCIIRI (6i)

Async (comodulation measure between channel A and B)




440 GAsync (global comodulation measure between channel A and B) CCTIRI (6i)

GAsync (global comodulation measure between channel A and B)




450 AminusB (A channel relationship to B channel) CCTRI (6i)

AminusB (A channel relationship to B channel)




451 BminusA (B channel relationship to A channel) CCTRI (6i)

BminusA (B channel relationship to A channel)




452 PDelta (B channel relationship to A channel) CCTIRI (6i)

PDelta (B channel relationship to A channel)




453 RatioA (Ratio of 2 streams from 1 channel) CCIIRI (6i)

RatioA (Ratio of 2 streams from 1 channel)




RatioB (Ratio of 2 streams from 1 channel)




455 RatioAB (Ratio of A to B in reward band) CCIIRI (6i)

RatioAB (Ratio of A to B in reward band)




456 RatioBA (Ratio of B to A in reward band) CCIIRI (6i)

RatioBA (Ratio of B to A in reward band)




458 D/S-Ratio (Ratio of A-B to A+B) CCIIRI (6i)

D/S-Ratio (Ratio of A-B to A+B)




Unity (1-Ratio of A-B to A+B)




480 ZCompAB (ZComposite) CCIIRI (6i)

ZCompAB (ZComposite)

EEGer4 Technical Manual




SingleA (one channel of data input)




603 SingleB (one channel of data input) CCIRRI (6r)

SingleB (one channel of data input)




612 Sum (sum of two channels of data input) CCIRRI (6r)

Sum (sum of two channels of data input)




613 Differ (channel A minus channel B) CCIRRI (6r)

Differ (channel A minus channel B)




634 PsyncABAB (comodulation measure AB twice) CCTRRI (6r)

PsyncABAB (comodulation measure AB twice)





668 MonitorA (1-ch + 12 monitors) CCMMMMMMMMMM (14a)

MonitorA (1-ch + 12 monitors)




MonitorAB (2-ch + 12 monitors)






ZCompAB (ZComposite)




700 ZAsymm (Zscore amplitude asymmetry) CCIRI (5)


ZAsymm (Zscore amplitude asymmetry)




ZCohere (Zscore coherence)




ZPhase (Zscore phase)







ZAbsPwrB (Zscore Abs Amp B)




ZRelPwrA (Zscore Rel Pwr A)






707 ZPRatioA (Zscore Power ratio A) CCIRI (5)

ZPRatioA (Zscore Power ratio A)




ZPRatioB (Zscore Power ratio B)




710 ZAsymm (Zscore amplitude asymmetry) CCIIRI (6i)

ZAsymm (Zscore amplitude asymmetry)



711 ZCohere (Zscore coherence) CCIIRI (6i)


ZCohere (Zscore coherence)



712 ZPhase (Zscore phase) CCIIRI (6i)


ZPhase (Zscore phase)

EEGer4 Technical Manual




ZAbsPwrA (Zscore Abs Amp A)


EEGer4 Technical Manual




EEGer4 Technical Manual




ZRelPwrA (Zscore Rel Pwr A)




EEGer4 Technical Manual



ZPRatioA (Zscore Power ratio A)




ZPRatioB (Zscore Power ratio B)



- Input
- stream
- Display
- Inhibit
- Allow
- Reward
- Game


1100 SingleA (one channel of data input) CCRRI (5)

SingleA (one channel of data input)




1101 SingleB (one channel of data input) CCRRRI (5)

SingleB (one channel of data input)



1110 Sum (sum of two channels of data input) CCRRI (5)


Sum (sum of two channels of data input)



1111 Differ (channel A minus channel B) CCRRI (5)


Differ (channel A minus channel B)

EEGer4 Technical Manual




SingleA (one channel of data input)

EEGer4 Technical Manual




1201 SingleB (one channel of data input) CCRRIIIRI (8)

SingleB (one channel of data input)



1210 Sum (sum of two channels of data input) CCRRRIRI (8)


Sum (sum of two channels of data input)



1211 Differ (channel A minus channel B) CCRRRIRI (8)


Differ (channel A minus channel B)

EEGer4 Technical Manual




1300 SingleA (one channel of data input) CCRRII (6i)

SingleA (one channel of data input)




SingleB (one channel of data input)




1310 Sum (sum of two channels of data input) CCRRII (6i)

Sum (sum of two channels of data input)




Differ (channel A minus channel B)




2000 SingleA (one channel of data input) CCMMR (5r)

SingleA (one channel of data input)




2001 SingleB (one channel of data input) CCMMR (5r)

SingleB (one channel of data input)




2010 Sum (sum of two channels of data input) CCMMR (5r)

Sum (sum of two channels of data input)




2011 Differ (channel A minus channel B) CCMMR (5r)

Differ (channel A minus channel B)




2020 Psync (synchrony measure between channel A and B) CCMMR (5r)

Psync (synchrony measure between channel A and B)




2030 Async (comodulation measure between channel A and B) CCMMR (5r)

Async (comodulation measure between channel A and B)




2040 GAsync (global comodulation measure between channel A and B) CCMMR (5r)

GAsync (global comodulation measure between channel A and B)




2050 AminusB (A channel relationship to B channel) CCMMR (5r)

AminusB (A channel relationship to B channel)




2051 BminusA (B channel relationship to A channel) CCMMR (5r)

BminusA (B channel relationship to A channel)




2052 PDelta (Peak Time Coherence) CCMMR (5r)

PDelta (Peak Time Coherence)




2060 ZCompAB (ZComposite) CCMMR (5r)


ZCompAB (ZComposite)





2100 ZAsymm (Zscore amplitude asymmetry) CCMMR (5r)

ZAsymm (Zscore amplitude asymmetry)











2104 ZAbsPwrB (Zscore Abs Amp B) CCMMR (5r)


ZAbsPwrB (Zscore Abs Amp B)






2106 ZRelPwrB (Zscore Rel Pwr B) CCMMR (5r)


ZRelPwrB (Zscore Rel Pwr B)





2108 ZPRatioB (Zscore Power ratio B) CCMMR (5r)


ZPRatioB (Zscore Power ratio B)



6000 QZcomp (4 channel zcomposite) CCCCCRRRR (8z)

QZcomp (4 channel zcomposite)


EEGer4 Technical Manual



6010 QABCD (4 channel amplitude) CCCCCRRRR (8z)

QABCD (4 channel amplitude)


EEGer4 Technical Manual



6020 QPSAvg (4 channel Average Psyncs) CCCCRRRR (8z)

QPSAvg (4 channel Average Psyncs)


EEGer4 Technical Manual



6100 DZcomp (4 channel zcomposite) CCCCCIRRI (8p)


DZcomp (4 channel zcomposite)

EEGer4 Technical Manual



6110 DZcompBPAB (4 channel zcompositeAB + single A) CCCCCIRRI (8p)


DZcompBPAB (4 channel zcompositeAB + single A)



6114 DZcompBPCD (4 channel zcompositeCD + single A) CCCCCIRRI (8p)


DZcompBPCD (4 channel zcompositeCD + single A)

EEGer4 Technical Manual




DSsingleA (4 channel single A)

EEGer4 Technical Manual




DSingleB (4 channel single B)

EEGer4 Technical Manual



DSingleC (4 channel single C)


EEGer4 Technical Manual



6150 DSingleD (4 channel single D) CCCCCIRRI (8p)

DSingleD (4 channel single D)


EEGer4 Technical Manual




6160 DPsyncABCD (4 channel two psync) CCCCCIRRI (8p)

DPsyncABCD (4 channel two psync)


EEGer4 Technical Manual



EEGer4 Technical Manual




EEGer4 Technical Manual



6182 DSumABSumCD (4 channel SumAB SumCD) CCCCCIRRI (8p)

DSumABSumCD (4 channel SumAB SumCD)


EEGer4 Technical Manual



6184 DSumABDiffCD (4 channel SumAB DiffCD) CCCCCIRRI (8p)

DSumABDiffCD (4 channel SumAB DiffCD)


EEGer4 Technical Manual



6186 DDiffABSumCD (4 channel DiffAB SumCD) CCCCCIRRI (8p)

DDiffABSumCD (4 channel DiffAB SumCD)


EEGer4 Technical Manual



6188 DDifffABDiffCD (4 channel DiffAB DiffCD) CCCCIRRI (8p)


DDifffABDiffCD (4 channel DiffAB DiffCD)

EEGer4 Technical Manual




DSingleAB (4 channel single A B)

EEGer4 Technical Manual




DSingleCD (4 channel single C D)

EEGer4 Technical Manual




EEGer4 Technical Manual




DSyncBPCD (4 channel PsyncCD+BP A)

EEGer4 Technical Manual




DPsync (4 channel PsyncAB+PsyncCD)

EEGer4 Technical Manual




DAsyncBPAB (4 channel AsyncAB+BP C)

EEGer4 Technical Manual



DSsyncBPCD (4 channel AsyncCD+BP A)


EEGer4 Technical Manual




6270 DAsync (4 channel AsyncAB+AsyncCD) CCCCIIRIIRI (10c)

DAsync (4 channel AsyncAB+AsyncCD)

EEGer4 Technical Manual




EEGer4 Technical Manual



DGAsyncBPCD (4 channel GAsyncCD+BP A)


EEGer4 Technical Manual




6300 DGAsync (4 channel GAsyncAB+GAsyncCD) CCCCIIRIIRI (10c)

DGAsync (4 channel GAsyncAB+GAsyncCD)

EEGer4 Technical Manual




EEGer4 Technical Manual



DPdeltaBPCD (4 channel PsyncCD+BP A)


EEGer4 Technical Manual



6330 DPdelta (4 channel PsyncAB+PsyncCD) CCCCIRIIIRI (10c)


DPdelta (4 channel PsyncAB+PsyncCD)

EEGer4 Technical Manual




DMinusBPAB (4 channel MinusAB+BP C)

EEGer4 Technical Manual



DMinusBPCD (4 channel MinusCD+BP A)


EEGer4 Technical Manual



6360 DMinus (4 channel MinusAB+MinusCD) CCCCCIRIIIRI (10c)


DMinus (4 channel MinusAB+MinusCD)

EEGer4 Technical Manual




DZcompBPAB (4 channel ZcompAB+BP C)

EEGer4 Technical Manual




EEGer4 Technical Manual




DZcomp (4 channel ZcompAB+ZcompCD)

EEGer4 Technical Manual




DRatioBPAB (4 channel RatioAB+BP C)


EEGer4 Technical Manual




DRatioBPCD (4 channel RatioCD+BP A)



EEGer4 Technical Manual




EEGer4 Technical Manual




DDifSumBPCD (4 channel Diff-SumCD+BP A)

EEGer4 Technical Manual




EEGer4 Technical Manual




DUnityBPAB (4 channel UnityAB+BP C)

EEGer4 Technical Manual




DUnityBPCD (4 channel UnityCD+BP A)

EEGer4 Technical Manual



DUnity (4 channel UnityAB+UnityCD)


EEGer4 Technical Manual



6600 QPSAvgA (4 channel PSync Averaged + 3 Inhibits) CCCCIIRI (8q)


QPSAvgA (4 channel PSync Averaged + 3 Inhibits)

EEGer4 Technical Manual




QSingleA (4 channel BP-A)

EEGer4 Technical Manual




QSingleB (4 channel BP-B)

EEGer4 Technical Manual




EEGer4 Technical Manual




QSingleD (4 channel BP-D)


EEGer4 Technical Manual




EEGer4 Technical Manual



EEGer4 Technical Manual




EEGer4 Technical Manual



6840 QQSsingleB (4 channel Monitor + Reward+Inhibits) CCCCCMMMMIRI (12)


QQSingleB (4 channel Monitor + Reward+Inhibits)

EEGer4 Technical Manual




QQSingleB (4 channel Monitor + Reward+Inhibits)

EEGer4 Technical Manual



QQSingleC (4 channel Monitor + Reward+Inhibits)

EEGer4 Technical Manual



QQSingleD (4 channel Monitor + Reward+Inhibits)

Appendix D: Data formats

Common Definitions:

```
/* common definitions */
/*
B4.0.0
030825.1126 hpl added definitions for devices
030905.1545 hpl defined c2mini
4.0.3g
040329.1650 hpl changed defs for multitude of J&J devices
4.0.99
040622.1406 hpl game defs
040830.1115 hpl reward mode group def
041026.1520 hpl other rsf threads (like tactile)
050422.1650 hpl added more device codes
050427.1420 hpl added display band type code
081104.1410 hpl Added FR_BAND
091004.1708 mjb renamed FF_DAK1 FF_UNITY
100824.1345 hpl added definition for eye status change
110209.1800 hpl pragma pack
110228.1700 hpl revisions for 4ch
*/
#ifndef __eegerh__
#define __eegerh__

#define TIMELOOKER 1

#pragma pack(push,1)

#define STATE_REMAP 0      // indeterminate
#define STATE_SETUP 1      // no drawing (not running)
#define STATE_SYNC 2       //
#define STATE_QUIT 9       // end program

#define STATE_DRAWABLE 10  // beginning of drawable states
#define STATE_INIT 10      // reinit state (starting a 'run')
#define STATE_FREEZE 11
#define STATE_PAUSE 12
#define STATE_RUN 20
#define STATE_REST 21
```

EEGer4 Technical Manual

```
// these are the queue numbers for the termination task queues
#define HSD_MAIN      1
#define HSD_THREAD    2
#define RSF_MAIN      3
#define RSF_THREAD    4
#define GAME_THREAD   5
#define RSF_OTHER     6

#define DISP_SETUP      0
#define DISP_HSD       1
#define DISP_LONGT     2
#define DISP_SPECT     3
#define DISP_ZSCORE    4
#define DISP_PERIODS   5

#define DISP_DEBUG     8
#define DISP_HELP      9

#define DISP_LSD      10 // there are up to 9 lsds
#define DISP_LSD1     11

#define DISP_EDITOR    19

#define DISP_GAME     20

#define SCRNMODE_SINGLE 2
#define SCRNMODE_DUAL   1
#define SCRNMODE_2COMP  0

#define RAWD 1          // raw data
#define PBACK 2         // playback
#define SIGGEN 3        // signal generator

#define DEV_PROSB    1
#define DEV_ADC      2
#define DEV_JJ       3
#define DEV_LOCKET   4
#define DEV_TTUSB     5
#define DEV_PET      6
#define DEV_NP24     7
#define DEV_BMSTR    8
#define DEV_QPET     9
#define DEV_A400      10
#define DEV_ATLANTIS 11

#define MAXUSEDSTREAMS 16 // most streams used for display (upper values used internally)
#define MAXPERIPHES  4 // maximum number of peripheral channels
```

EEGer4 Technical Manual

```
#define MAXPERIPHSALLOWED 3
#define MAXLOWPASS 4      // maximum number of lowpass channelsstreams
#define FLT_R1    MAXUSEDSTREAMS // offset to backing store for reward 1
#define FLT_R2    (FLT_R1+MAXUSEDSTREAMS) //offset to backing store for reward 2
#define FLT_R3    (FLT_R2+MAXUSEDSTREAMS) //offset to backing store for reward 3
#define FLT_R4    (FLT_R3+MAXUSEDSTREAMS) //offset to backing store for reward 4
#define FLT_LIVE   (FLT_R4+MAXUSEDSTREAMS) // offset to live data for sham runs
#define FLT_NOTCH  (FLT_LIVE+MAXLOWPASS+MAXPERIPHS) // offset to notch filter coefficients
#define FLT_NOLO   (FLT_NOTCH+MAXPERIPHS) // offset for non-lowpass, non-dc-corrected raw values (for SCP and VLF)
#define FLT_PERIPH (FLT_NOLO+MAXLOWPASS) // offset for peripheral state values which may not be used...
#define FLT_SIZER   (FLT_PERIPH+MAXPERIPHS) // allocation counter for states

#define MAXPERIODS 128      // maximum number of periods in a session
#define MAXREWARDS 4        // maximum number of reward streams
#define BACKINGSETS 4
#define MAXDEVICES 2        // maximum number of input devices
#define MAXDEVCHANNELS 4    // maximum number of highspeed channels per device
#define MAXSTRANDS 8         // strands in EGS games
#define MAXREWARDSTRANDS 4   // maximum number of reward strands
// message ids
#define EEGMBASE 0x1000

#define EMSG_RSFTIME WM_APP+EEGMBASE+ 0    // hs timer message
#define EMSG_RSFWDOG WM_APP+EEGMBASE+ 1    // sent to OC to show alive

#define EMSG_HSDTIME WM_APP+EEGMBASE+100    // time to run
#define EMSG_HSDWDOG WM_APP+EEGMBASE+101    // sent to OC to show alive

#define EMSG_OCWDOG WM_APP+EEGMBASE+201    // sent by OC to RSF to show alive
#define EMSG_OCACT  WM_APP+EEGMBASE+202    // operator control message to RSF

// the GROUP codes

#define GROUP_DC    0    // Display/control settings
#define GROUP_SCALE 1    //Scale settings
#define GROUP_MODE  2    //Mode changes (up/down, model, etc)
#define GROUP_RWD   3    //reward mode changes
#define GROUP_SITE  4    //Site location change
#define GROUP_STAGE 5    // Stage setting
#define GROUP_THRESH 6   //Thresh settings
#define GROUP_LOFR  7    // low freq or center freq
#define GROUP_HIFR  8    // high freq or width
#define GROUP_AGOAL 9    // autoset (1 bit), % (7 bits), time (byte), min (byte), max (byte)
#define GROUP_MISC  14    // some kind of setting to record (like integration time)
```

EEGer4 Technical Manual

```
#define GROUP_EVENT 15 // some kind of event

// the ACTION codes
#define ACTION_SET 00 // just set the value
#define ACTION_KEY 01 // keycode
#define ACTION_MISC 02 // misc meaning

#define ACTION_REWARD 15 // reward granted
#define ACTION_FBACK 16 // begin feedback
#define ACTION_PAUSE 17
#define ACTION_USER 18 // user event (f8)
#define ACTION_USERM 19 // user event with message
#define ACTION_QUIT 20 // quit session

#define MAKE_ACTIONCODE(group,action,stream) (((group & 0x0f) << 12) | ((action & 0xff) << 4) | (stream & 0x0f))
#define ACTIONCODE_GROUP(ac) ((ac >> 12) & 0x0f)
#define ACTIONCODE_ACTION(ac) ((ac >> 4) & 0xff)
#define ACTIONCODE_STREAM(ac) (ac & 0x0f)

// the stage kind codes
#define ST_SETUP 0
#define ST_RUN 1
#define ST_PAUSE 2
#define ST_EXIT 3
#define ST_PERIPH 4
#define ST_BASELINE 5

// the stage sequence codes
#define SEQ_INIT 0
#define SEQ_PAUSE 1
#define SEQ_RUN 2
#define SEQ_REST 3
#define SEQ_STAGEP 4
#define SEQ_CHANGE 5

// message block structure
typedef struct _mymsgblock_
{
    short code;      // why there is a block
    short len;       // number bytes INCLUDING block headers
    long filler;     // forcing to 16
    union {
        double d;
        long l[2];
        short s[4];
        char c[8];
    } u;
}
```

EEGer4 Technical Manual

```
 } MYMSGBLOCK;

#define MB_KEY      1      // means a keycode in s[0];
#define MB_ACTION   2      // means an action code in s[0];
#define MB_KEYBD    3      // WM_CHAR keys for general typing

#define MBMOD_SHIFT 1      // shift key
#define MBMOD_CTRL  2      // ctrl key
#define MBMOD_ALT   4      // alt key

typedef struct _fcmd_ {

    unsigned char funct;           //Function code     8
    unsigned char flags;           //flag bits        8
    // rrxxgemd
    //          |---- 'Display' this stream
    //          |---- multiple inputs
    //          |----- more than up/down rewards allowed
    //          |----- outputs to game
    //          |
    //          |
    //          ||----- backing set (0-3)selects FLT_R1->FLT_R4
    //
    unsigned char grpstream;
    // ggssssss
    //      ||||||----output stream
    //      ||-----reward group
    unsigned char procmode;         // code for proc mode
    char numinp;                  // number of inputs 8
    char inch[3];                 // input array      24      up to 6 inputs packed as seq of bytes
    char numinhib;                // number of inhibits 8
    char inhib[3];                // array            24      up to 6 inhibits packed as seq of bytes
} FCMD;
#define FCMD_FLAG_DISPLAY 1
#define FCMD_FLAG_MULTINP 2
#define FCMD_FLAG_EXTRA   4
#define FCMD_FLAG_GAME    8

// stream usage codes
// These orders must match parblock.py !!!!
#define SUSE_RAW      1
#define SUSE_INHIB    2
#define SUSE_REWARD   3
#define SUSE_MONITOR  4
#define SUSE_DISPLAY  5

#define MAX_FILTER_MODES 32           // in a single session
```

EEGer4 Technical Manual

```
#define MAX_FM_FUNCTIONS      32          // most filter functions per layout

// filter function names
// These names/orders must match parblock.py !!!!
#define FF_END              0          //END of function list
#define FF_LOWPASS           1          //lowpass
#define FF_BANDPASS           2          //bandpass
#define FF_DIFFER              3          //difference
#define FF_SYNCH              4          //synch reward
#define FF_COMOD              5          //comod reward
#define FF_GLCOM              6          // global comod
#define FF_MINUS              7          // XminusY
#define FF_SMRREW             8          //single SMR reward
#define FF_ATREW              9          //single AT reward
#define FF_COHERE             10         // coherence mode
#define FF_MULTIREW           11         // multi-reward
#define FF_RATIO               12         // ratio
#define FF_DAKCOH              13
#define FF_DAKPHASE            14
#define FF_ZSCOREAA             15
#define FF_ZSCORECO             16
#define FF_ZSCOREPH             17
#define FF_ZSCOREAPA            18
#define FF_ZSCOREAPB            19
#define FF_ZSCORERPA            20
#define FF_ZSCORERPB            21
#define FF_ZSCOREPRA            22
#define FF_ZSCOREPRB            23
#define FF_DC                  24
#define FF_UNITY                25
#define FF_DAK2                 26
#define FF_DAK3                 27
#define FF_DAK4                 28
#define FF_DIFFSUM              29          // added for Hirschberg study
#define FF_ZCOMPOSITE           30

#define FR_UP                  0          // normal "UP" reward
#define FR_DOWN                1          // down reward
#define FR_CENTER               2
#define FR_ELEVATE              3
#define FR_DEVUP                4
#define FR_DEVDOWN              5
#define FR_BAND                 6
#define FR_ZCOMP                 7
```

EEGer4 Technical Manual

```
#define GMISC_INTEGTIME          1      // filter integration time
#define GMISC_TIMEONSTATE         2
#define GMISC_PERCENTONSTATE      3
#define GMISC_MINREWARD_TIME      4
#define GMISC_EYESSTATUS          5

#define GMISC_ZCPARAMETER 16
#define GMISC_ZCHAN        17
#define GMISC_ZCBAND        18
#define GMISC_ZCTHRSH       19
#define GMISC_ZCPERCENT     20
#define GMISC_ZCENABLE       21
#define GMISC_ZCDISABLE      22

#pragma pack(pop)
#endif
```

Raw Data File Format

```
/* file format header */

/*
RAW data files
The raw filename is constructed from the client code (<32 characters, no embedded spaces
or special characters) and a numerical code. As with the summary filename, a 6 digit
code will be APPENDED to the client code. The 6 digits are the 4-digit partial Julian
daycode and a uniqueness-guaranteeing sequence code.
The filename structure will be CCffffqq.RAW where CC is the 1- to 31-character client
id code, ffff is the part of the Julian date, and qq is the sequence number.
```

Raw File Format

Basically, a 'raw' file contains 'raw' data acquired during a session along with operator actions time-correlated with the data. Up to 16 channels of data are provided for, each with scaling, format definitions, etc. Channel data is stored end-to-end, not interleaved. There may be multiple blocks of channel data for each channel in the file. Multiple blocks will be time-sequential in the file although the block positions imply no time relationships between channels. By this, I mean that a low speed block may contain timed samples over a longer time period than higher-rate blocks before and after it sequentially in the file. The actual I/O may/will be

EEGer4 Technical Manual

double-buffered to handle the 'long' time it takes to write data blocks to a file.

Header format	Content	Symbol	Type	Length	Notes
File type code		char		4	RAWD
Format code		short int	2 bytes		Format code is major.minor , major* 100 format 0.1 == 1, 1.3 == 103
Datecode		long int	4		Abs date
Timecode		long int	4		Seconds since midnight of FB start time
Client code		char	32		client id code
Client name		char	64		client name
Machine ID		long int	4		dongle code for now
Sample Clock Rate		short int	2		e.g. 160, 256, 512
number of channels		short int	2		1 to 16
format string		char	16		
offset of controls					
block		long int	4		byte offset
Max datablock size		long int	4		size of largest data block must allocate (1+ #chans)*2 of these

Chan0Block 16 channel block

Chan1Block

Chan2Block

Chan3Block

Chan4Block

Chan5Block

Chan6Block

Chan7Block

Chan8Block

Chan9Block

Chan10Block

Chan11Block

Chan12Block

Chan13Block

Chan14Block

Chan15Block

Format of each channel block

Content	Symbol	Type	Length	Notes
Type code		unsigned	2	drives format of data
substrate		unsigned	2	currently only 1 or 8 for procomps
SequenceCode		long	4	unique sequence code for stream
ID		long	4	procomp serial number or ??
Scale factor		double	8	lsb scale factor of data
offset of 1st data block		long	4	byte offset

EEGer4 Technical Manual

Format of a data block

Content	Symbol	Type	Length	Notes
Length of block		long int	4	length in bytes of block incl header
Type code			2	0=uncompressed
time stamp of 1st data		long int	4	real-time cycle count
offset of next block		long int	4	byte offset of next block
data			2	
----			2	

Format of controls block

Content	Symbol	Type	Length	Notes
Length of block		long int	4	length in bytes of block incl header
offset of next control blk		long	4	byte offset of next block

Format of control sub-blocks (each corresponding to one action!)

length		2 bytes	sub-block length including this header
action code		2 bytes	LOTS of codes
		Code definitions use lots of major grouping codes for easier categorization.	
		Lower 4 bits of every code reflects (possible) channel.	
time stamp		4 bytes	real-time cycle count
data		double	
		char [8]	
		long[2]	8 bytes
this may be the beginning of a variable length text field also			
(for event reasons and such)			

Action coding:

Multiple actions with same time stamp imply a major action (such as beginning a session using the previous settings!!).

MMMAAAAAAAACCCC

MMMM is the major grouping
00 Display/control settings
01 Scale
02 Mode changes (up/down, model, etc)
03
04 Site location change
05 User event
06 Threshold setting

EEGer4 Technical Manual

```
07 Low freq setting
08 high freq setting

AAAAAAA is the action code
    tbd - defined in eeger.h
CCCC is the channel for channel-related action codes
060214.1135 hpl reverted to 104 version temporarily
070822.1040 hpl made some shorts into unsigned short
100621.1617 hpl added subset code to header and biological sex
110228.1500 hpl 110 revised header layout and made room for more peripheral data
*/
#ifndef __rawformat_
#define __rawformat_
#pragma pack(push,1)

#define RAWVERSION 110      // odd versions are the 'compressed' version
                        // of the next lower even file format
                        // this makes it easy to discriminate/convert

#define MAX_CHUNKS 64        // maximum number of data chunks in a file
#define MAX_ACTION_SIZE 4096 // max size of an action chunk

// there are 16 of these in the header block
typedef struct _RAWCHD_
{
#if RAWVERSION >= 110
    unsigned char typecode;      //my code for device type for eeg data (DEV_ codes
    unsigned char channel;       // channel number for data A=0, B=1, C=2, D=3, whatever for peripherals, etc.
#else
    unsigned short typecode;     //which device produced data
#endif
    unsigned char substrate;     //currently only 1 or 8 for procomps
                                // Note: substrate == 1 means raw EEG data
                                //          substrate != 1 ==> peripheral data
    unsigned char datatype;      // kind of data
                                // 0 == EEG, 1= undifferentiated peripheral data, 2->255 ???
    long ID;                   //procomp serial number or peripheral channel
    double scalefactor;         //lsb scale factor of data
                                // peripheral data has a scalefactor of ????

    long dboffset;              //byte offset to first data block
} RAWCHD;

// this is the actual file header
typedef struct _RAWHD_
{
```

EEGer4 Technical Manual

```
char filetype[4];           //RAWD
short formatcode;           //Format code is major.minor ,major* 100
                           //format 0.1 == 1, 1.3 == 103
long datecode;              //Abs date
long timecode;              //Seconds since midnight of FB start time
char clientcode[32];         //client id code
char clientname[64];         //client name
long machineID;             //dongle code for now
short clockrate;             //e.g. 160, 256, 512
short numberofchannels;      //1 to 16
char formatstring[16];        // contains format string at recording time
long SequenceCode[MAX_FILTER_MODES]; // sequence codes for session !!!
long cboffset;               //byte offset of controls block
long maxdatablocksize;       //size of largest data block
                             //must allocate (1+ #chans)*2 of these

#if RAWVERSION < 110
    RAWCHD chd[15];          //channel head blocks
    char lcode;                // 0 for unk, 1 for therap, 2 for remote
    char subset;                // which subset of current structure this is
    char bsex;                  // either 0 (unk), 'M','F'
    char filler[17];           // make up for stealing one chd
#endif
#if RAWVERSION >= 106
    long birthdate;            // birthdate in "proleptic Gregorian" == "absolute date" format 1= 01/01/0001
#endif
#if RAWVERSION >= 108
    char xguid[32];            // 'original' writer of file
#endif
#else
    char lcode;                // 0 for unk, 1 for therap, 2 for remote, 3 for sham
    char subset;                // which subset of current structure this is
    char bsex;                  // either 0 (unk), 'M','F'
    char filler[17];
    long birthdate;            // birthdate in "proleptic Gregorian" == "absolute date" format 1= 01/01/0001
    char xguid[32];            // 'original' writer of file
    long chdoffset;             // byte offset to first (of the sequential) channel block
    long summaryoffset;          // byte offset to included summary file if there is one appended!!!!
    RAWCHD chd[16];             //channel head blocks
#endif
} RAWHEAD;

typedef struct _OLDRAWHD_
{
    char filetype[4];           //RAWD
short formatcode;           //Format code is major.minor ,major* 100
                           //format 0.1 == 1, 1.3 == 103
long datecode;              //Abs date
long timecode;              //Seconds since midnight of FB start time
char clientcode[32];         //client id code
```

EEGer4 Technical Manual

```
char clientname[64];           //client name
long machineID;                //dongle code for now
short clockrate;                //e.g. 160, 256, 512
short numberofchannels;         //1 to 16
char formatstring[16];          // contains format string at recording time
long SequenceCode[MAX_FILTER_MODES]; // sequence codes for session !!!
long cboffset;                 //byte offset of controls block
long maxdatablocksize;         //size of largest data block
//must allocate (1+ #chans)*2 of these
RAWCHD chd[15];                //channel head blocks
char lcode;                     // 0 for unk, 1 for therap, 2 for remote
char subset;                    // which subset of current structure this is
char bsex;                      // either 0 (unk), 'M','F'
char filler[17];                // make up for stealing one chd
#if RAWVERSION >= 106
    long birthdate; // birthdate in "proleptic Gregorian" == "absolute date" format 1= 01/01/0001
#endif
#if RAWVERSION >= 108
    char xguid[32];           // 'original' writer of file
#endif
#endif

} OLDRAWHEAD;

// there is one of these at the head of each chunk of raw data
typedef struct _RAWDB_
{
    long length;                  //length in bytes of block incl header
    unsigned char typecode;        //device it came from
    unsigned char datacode;        // how to decode special data
    long timestamp;                //real-time cycle count of first sample
    long nextoffset;               //byte offset of next block
    //short data[1];                // data sample
    // ...
    more samples
} RAWDB;

// there is one of these for each 'action' subblock
typedef struct _RAWSUB_
{
    unsigned short length;         //sub-block length including this header
    unsigned short actioncode;      //LOTS of codes
    long timestamp;                //real-time cycle count of this action
    union {
        double data1;              //this may be the beginning of a variable length text field also (for event reasons and such)
        char text[8];
        long ld[2];
        float fdata[2];
        short twobytes[4];
    } d;
}
```

EEGer4 Technical Manual

```
    } RAWSUB;

    // there is one of these for each chunk of actions
    typedef struct _RAWCB_
    {
        long length;                      //length in bytes of block incl header
        long nextoffset;                  // offset in bytes of next control blk
        //RAWSUB blocks[1];                // sub-blocks.
    } RAWCB;

#pragma pack(pop)
#endif
```

Summary File Format

```
/*
Summary files consist of a header and a series of data blocks containing the once per second
average values, the threshold settings, and a composite reward flag bit for each filtered
channel. Additional information is saved (at the slow 1 Hz rate) for later
analysis/understanding
```

```
Part of the summary data filename is constructed from the 4 lower digits of the Julian date.
Julian dates in the range 1995-2023 are all greater than 245000 and less than 255000.
Additionally, (to handle multiple and aborted sessions in the same day),
a 2-digit sequence number will be used to ensure that there are no filename conflicts.
The filename structure will be Sjjjjqq.SUM where jjjjj is the aforementioned part of the
Julian date and qq is the sequence number. Internally, the actual client identifiers are
stored in the file.
```

Summary file data format

This file consists of a header and a series of data blocks containing the
once per second average values, the threshold settings, and a composite
reward flag bit for each filtered channel.
Additional information is saved (at the slow 1 Hz rate) for later
analysis/understanding.

For data blocks				
Content	Symbol	Type	Length	Notes

EEGer4 Technical Manual

```
0th average value      short int    2      microvolts * 100; negative means reward for at least one sample in 1 second interval
        Note that maximum average microvolt value
        is 327.67 microvolts (times 100!!)
..
nth  average value      n is number of traces -1
0th average value      for the next second!!!
-----
For threshold blocks
Content      Symbol  Type      Length  Notes
0th threshold      unsigned   2      microvolts * 100

..
nth  threshold      n is number of traces -1

For frequency blocks
Content      Symbol  Type      Length  Notes
0th low freq      unsigned   2      low freq in 1000ths of Hz; 4250 = 4.25 Hz
0th high freq     unsigned   2      in 1000ths of Hz
-----
nth low freq      n is number traces - 1
nth high freq

Summary block
Content      Symbol  Type      Length  Notes
0th average value      coded    2      microvolts * 100
        Note that maximum average microvolt value
        is 327.67 microvolts (times 100!!)
..
0th percentage      in % times 100
-----
nth percentage
nth  average value      n is number of traces -1

Scale block
Content      Symbol  Type      Length  Notes
0th scale      unsigned   2      microvolts * 100

..
nth  scale      n is number of traces -1
```

EEGer4 Technical Manual

```
030326.1151 hpl added mark data but kept version 106
040127.1152 hpl changed to unsigned short length and version 108
4.1.xx
040913.1130 hpl version 110 - added overall reward percent to period data
050420.1115 hpl version 112 - scab on total reward % to percentage values (fake extra trace data)
060120.1400 hpl version 114 - add SUM_PERIPH to data, go to long header lengths
060214.1135 hpl reverted to 112 structures temporarily
070129.1437 hpl 116 adds zscore data
110209.1800 hpl pragma pack
110228.1510 hpl 118 revised structure to remove xdata and handle both 2- and 4-channel zscore
*/
#ifndef __SUMFORMAT_H__
#define __SUMFORMAT_H__
#include "eeger.h"
#pragma pack(push,1)

#define SUMVERSION 118

// file header
typedef struct _SUMHD_
{
    char filetype[4];           // SUMD
    short formatcode;          //Format code is major.minor ,major* 100
                               //format 0.1 == 1, 1.3 == 103
    char protocolcode[4];       //Like SMR AT EXP with trailing null bytes
    char numberoftraces;       //1 to 16 filtered traces, 1-16 channels
    char numberlowpass;        // number lowpass channels
    long datecode;              //Abs date
    long timecode;              //Seconds since midnight local time
    long SequenceCode[MAX_FILTER_MODES];// what kind of sequence codes were possible
    char formatstring[MAXUSEDSTREAMS];    // contains format string at recording time
    char gameid[30];
    char lcode;                  // 1 for therapist, 2 for remote
    char spare;
    char clientcode[32];
    char clientname[64];
    // following new in 118
    char numberperiph;          //number of peripheral channels
    float multiplier[MAXUSEDSTREAMS+MAXPERIPH]; // scaling factor for data - mostly peripherals!!
} SUMHD;

typedef struct _SUMBLKHD_
{
#if SUMVERSION >= 114
    unsigned long length;        //length of this block in bytes
#else
    unsigned short length;

```

EEGer4 Technical Manual

```
#endif
    short formatcode;           // SUM_ codes
    unsigned short timestamp;   // seconds relative to base timecode of 1st data
} SUMBLKHD;

typedef struct _SUMFR_
{
    SUMBLKHD blkhd;
    struct {
        unsigned short low; // low freq in 1000ths of Hz; 4250 = 4.25 Hz
        unsigned short high; // in 1000ths of Hz
    } freq[1];

    // nth low freq          n is number traces - 1
    // nth high freq
} SUMFREQ;

typedef struct _SUMCD_
{
    SUMBLKHD blkhd;
    char modename[16];          //kind of feedback mode
    struct {
        char channel;          //channel for input
        char sites[11];         //string of sites
    } site[1];
    //Chan code             pairs of entries for all USED channels
    //Site code
} SUMCD;

//missing data points denoted by value of 0xffff
typedef struct _SUMDTA_
{
    SUMBLKHD blkhd;
    union {
        short value;           //microvolts * 100; negative means reward for at least one sample in 1 second interval
                                // Note that maximum average microvolt value
                                // is 327.67 microvolts (times 100!!)
        unsigned short threshold;//microvolts * 100
        unsigned short average; //microvolts * 100 (long term average)
                                //Note that maximum average microvolt value is 327.67 microvolts (times 100!!)
        unsigned short percentage; // % times 100
        // I add an EXTRA value == total reward percentage on
        unsigned short scale;    // microvolts * 100
        unsigned short periph;   // something * 10
    } trace[1];
    // nth percentage
```

EEGer4 Technical Manual

```
//nth    average value          n is number of traces -1
} SUMDATA;

typedef struct _SUMPER_
{
    SUMBLKHD blkhd;
    unsigned short period;           // period number
    unsigned short seconds;          // delta seconds in period
    unsigned short score;            // delta score for period
    unsigned short overall_reward;   // overall reward % new in 110
    char gamename[64];              // rightmost 64 chars of game 'name'/path

// begin junk
// unused since at least 4.2.0
    struct {
        unsigned short beginampl;    // 1st second ampl * 100
        unsigned short aveampl;      // average amplitude * 100
        unsigned short endampl;       // end of period ampl * 100
        unsigned short lowfreq;       // band low freq * 1000
        unsigned short highfreq;      // band high freq * 1000
    } xbands[4];
// end junk
    struct {
        unsigned short average;       // ending long-term average * 100
        unsigned short percent;       // ending percent over threshold *100
        unsigned short threshold;     // ending threshold value * 100
        unsigned short lowfreq;       // ending low freq limit * 1000
        unsigned short highfreq;      // ending high freq limit * 1000
        char sitelist[12];           // ending site list
    } streamval[1];
} SUMPERIOD;

typedef struct _SUMMODE_
{
    SUMBLKHD blkhd;
    long seqcode;                  // sequence code now switching to
    char modename[16];             // kind of feedback mode
    unsigned char rewardcode[16];   // codes for reward directions
} SUMMODE;

typedef struct _SUMMARK_
{
    SUMBLKHD blkhd;
    unsigned short msglength;      // number of bytes in msg
    char msg[1];                  // where mesg starts
} SUMMARK;

typedef struct _SUMZSCORE_
{
```

EEGer4 Technical Manual

```
SUMBLKHD blkhd;
short int kindcode;                                // 0 means all 2ch terms are included;1 means only products, 2= all 4ch,3 = only 4ch products
union __kinddata__
{
    struct __chans4__
    {
        struct {
            short int zaa[6][10];      // microvolts*100
            short int zco[6][10];
            short int zap[6][10];
        } products[1];
        struct {
            short int zap[4][10];      // 1st 8 are chan A
            short int zrp[4][10];
            short int zpr[4][10];      // 1st 10 are chan A
        } terms[1];
    } chans4;
    struct __chans2__
    {
        struct {
            short int zaa[1][10];      // microvolts*100
            short int zco[1][10];
            short int zap[1][10];
        } products[1];
        struct {
            short int zap[2][10];      // 1st 8 are chan A
            short int zrp[2][10];
            short int zpr[2][10];      // 1st 10 are chan A
        } terms[1];
    } chans2;
} kinddata;
} SUMZSCORE;

//WARNING: storage.cpp and sumfile.py (and others!?!?)
// make assumptions about these numbers!!!!
#define SUM_DATA      0      //SUMDATA
#define SUM_THRESH     1      //SUMDATA
#define SUM_AVERAGE    2      //SUMDATA
#define SUM_PERCENT    3      //SUMDATA
#define SUM_FREQ       4      //SUMFREQ
#define SUM_SCALE      5      //SUMSCALE
#define SUM_SITE       6      //SUMCD
#define SUM_PERIOD     7      //SUMPER
#define SUM_MODE       8
#define SUM_MARK       9
#define SUM_PERIPH    10 //SUMDATA
#define SUM_ZSCORE     11
#define SUM_LAST       11      // just the highest number
```


EEGer4 Technical Manual

```
#pragma pack(pop)
#endif
```

Appendix E: Data Acquisition Methodology

This document contains data acquired using normal versions of EEGeri4. The data is generated using a combination of a software and hardware external to EEGeri4 itself as described in Appendix F.

EEGeri4 supports a special development mode where a specific run stage comment in a session plan sends a UDP message to the signal generator described in Appendix F. The message controls some EEGeri4 parameters and the frequency that the signal generator is to supply. Critical parameters other than the frequency are the settings for the lowpass filter and DC correction filters which can be turned on/off depending on the test to be run. The EEGeri4 setup is done in session planning where the alert name must be SENDUDP.WAV and the remark field contains the parameters.



Frequencies used in the test session plan are

0.0,0.1,0.2,0.3,0.5, Hz

1 to 40 Hz in 1 Hz steps

45.0,47.0,48.0,49.0,50.0,51.0,52.0,58.0,59.0,60.0,61.0,62.0,65.0,70.0 Hz

The run time chosen is long enough for the value smoothing to occur so that values at the end of each period represent the average value. Since there is no pause flag checked, run stages continue to run until the end of the session. In this fashion, a test run is semi-automated. Once started, the identical test sequence is performed. Data for each run is available as summary data outputs from the EEGeri4 review process.

This method was used to generate the filter characteristics used for the filter plots and the bandpass characteristics shown for each amplifier/encoder.

Appendix F: Signal Generator

Signal Generator

Technical Description

Version 1.12

Copyright © 2012 EEG Software LLC

Description

The signal generator tool provides a representative EEG signal for tests of signal acquisition devices and for end-to-end testing of the entire neurofeedback system.

Digital to Analog Converter (DAC)

The DAC used for the signal generator tool is a USB-3101FS 16-bit four channel DAC manufactured by Measurement Computing Corporation. It is supported by the MCC-provided InstaCAL input/output interface library (which must be installed). The device is a USB-2.0 controlled and powered device.




Figure 2:

Figure 1:


Screw terminal pin assignments

Terminal	Signal
0	AO 0
1	Common (COM)
2	AO 1
3	Common (COM)
4	AO 2
5	Common (COM)
6	AO 3
7	Common (COM)
8	NC (No connection)
9	Common (COM)


Attenuator

The attenuator device is an internally manufactured purely resistive voltage divider used to reduce the DAC outputs (0.5 volt) to ranges appropriate for an EEG-level (50 microvolt) sensing device.

Attenuator box.




Internals of typical attenuator box.




EEGer4 Technical Manual

Schematic of attenuator box.




Cabling from USB-3101FS to attenuator box.



EEGer4 Technical Manual

Wiring diagram for cable between USB-3101FS and attenuator box.



SigDriver

This software module contains the logic necessary to decode an input data file (in EEGer .RAW format) and send it to the DAC. It also provides the ability to generate simple sine wave signals and to provide a spike impulse for testing of EEGer software. The user interface for the module is implemented using FLTK, a user interface library. Source code for this tool is maintained in the revision control system. The module contains the following logical functions:

- mymain – builds the screen interface and calls the output routine
- runsome - the actual output routine
- findfiles – determine which files are available
- read_raw_file – reads/decodes a raw file
- read_first_event – reads/decodes first text event in a raw file (for titles)
- synth_raw_file – creates equivalent of raw data in cases of simple sine waves (manual) mode
- sample – the actual computation of the output voltage for a sample
- CheckListener – support routine for automatic bandpass testing which monitors for frequency changes commanded by EEGer in test mode.

Because this program runs under Microsoft Windows, the computation of the correct signal value takes into account the erratic timing under Microsoft Windows and handles it appropriately.

The format of the data listened for is a UDP message received containing one block of data in this format:

0xe1	0	Size low	Size high
Text			


Where “text” is an ASCII string in the following format:

freq:optn:optn:optn..

where only the first value (frequency value terminated by a colon) is used by SigDriver. This value is forced into the manual frequency setting for Sig1. At that time, all other values are set to zero except the amplitude for Sig1 (set to 40) and the mask for Sig1 (set to 3). Manual mode is forced on at that time.


EEGer4 Technical Manual

The main screen is this:




This reports that the testband file is being used as the signal source, there are two channels of data in the source, we are NOT forcing all the data channels to use one channel of data, and that the signal driver is running.

The dropdown box at the top (File to play) gives the following options (unless additional files are created using the FreqFileGen tool described later):




EEGer4 Technical Manual

Selection of the Manual checkbox gives this screen:




where the yellow-highlighted region contains the values to generate. Shown here in this example is a 10 Hz signal on channel 0 (A) and a 12 Hz signal on channel 1 (or B).

If the Dev checkbox is selected, the following screen shows the additional options for the impulse generator:



Frequency File Generator (*FreqFileGen*)

This module is a Python script which allows creation of complex source signal files (in EEGer format suitable for replay in EEGer). The screen looks like this:



There are 32 source points, each of which can have the following characteristics:

- frequency
- peak-peak voltage in microvolts
- starting time offset of signal (in seconds)
- amplitude variation frequency (rate of variation of amplitude)
- maximum amplitude of amplitude variation (amount of amplitude variation)
- jitter frequency (rate of variation of starting time)
- jitter max (maximum amount of time variation)
- Channel mask (bit mask for the output channels for the signal)
- Start (starting time of the signal)
- length(length of time signal is generated – 0 means forever)

In addition, each of the four possible channels can have a DC offset, a DC offset variation, and a starting time for the DC offset variation to test DC characteristics of AC and DC coupled acquisition devices.

EEGeri4 Technical Manual

Calibration process

The signal generator tool must be periodically calibrated. The calibration tool used is a Fluke 87-5 digital multimeter (DMM) or equivalent which has a valid NIST-traceable calibration certificate.

1. DISCONNECT signal block from USB-3101FS

Frequency/amplitude of signal generator:

2. The DMM is connected to the terminal block of the USB-3101FS pins 0 (signal) and 1 (ground).
3. The SigDriver is started and Manual mode is selected.
4. The following parameters are sequentially entered in the manual fields and results verified:

Frequency	Amplitude	Freq reading	Amplitude reading (true RMS)
2	10	2+- 0.3	0.10 +- .0 1
5	15	5 +- .3	0.15+- .01
10	20	10+- .3	0.20+- .01
12	25	12 +- .3	0.25+- .01
25	30	25+- .5	0.30+- .01
30	40	30 +- .5	0.40+- .01
50	60	50 +- .5	0.60+- .01

5. Repeats steps 2,3, and 4 for DMM connected to pins 2,4,6 making sure that the signal is routed to the correct channels (mask=15)

This verifies the generation of signals has requisite accuracy.

Attenuator box:

5. Using DMM, measure resistance between signal block pin 0 and signal block pin 1 while attenuator box is connected to cable but NOT to an acquisition device.
6. Verify resistance is 1,000,000 +- 20,000 ohms.
7. Measure resistance between signal and reference leads for channel A on the attenuator box.
8. Verify resistance is 100 +- 3 ohms
9. Repeat steps 5-6 for pins 2 and 1 on the signal block.
10. Repeat steps 7 and 8 for leads for channel B on the attenuator box.

Each attenuator box needs to be calibrated and a calibration sticker attached.

11. Fill in the calibration test log and sign where appropriate.

This procedure needs to be once per year or sooner if repairs are needed on any component or the driver software is changed.

EEGer4 Technical Manual

Initial calibration data

